

Liquid Culture Re-Amplification Protocol by Chong Park

For sublibrary amplification:

- 1) Mix: 10-50 ng of DNA + 100 uL DH5 α high efficiency (3×10^9 cfu/ug) (1 rxn per sublibrary)
- 2) Transform using the manufacturer's suggested protocol, making sure to perform the full 30min ice-incubation.
- 3) Recover in 1 mL total with SOC for 1 hour while shaking at 37C.

- 4a) Take 5ul of recovery. Make serial dilutions and plate them onto LB/Amp plate to calculate transformation efficiency.
- 4b) Add the rest of recovery to 500 mL LB+Carb. Grow O/N while shaking @ 37°C (16h)

- 5) Calculate transformation efficiency next day. If the efficiency is higher than 100~200 colonies per sgRNA construct in the library, harvest cells and purify the library.
- 6) Depending on your pellet weight use multiple Maxiprep columns (Sigma or Zymogen), a Megaprep (Qiagen or Sigma), or Gigaprep (Qiagen, Sigma, Zymogen).
- 7) Expect a yield ~2mg.

For pooled library

Use MegaX DH10B ($>3 \times 10^{10}$ cfu/ μ g) electrocompetent cells (Fisher cat# C640003) to achieve higher coverage (~1000 colonies/sgRNA).

- 1) Mix 100ng of pooled library with 50ul of MegaX cells.
- 2) Transform using the manufacturer's suggested protocol (2.0kV, 200 ohms, 25uF in 0.1cm cuvette). Do not ice incubate.
- 3) Recover in 1mL total with SOC for 1.5-2 hours while shaking at 37C
- 4a) Take 5ul of recovery. Make serial dilutions and plate with beads to calculate transformation efficiency.
- 4b) Add the rest of recovery to 500 mL LB+Carb. Grown O/N while shaking @ 37C (16h)

- 5) Calculate transformation efficiency next day. If the efficiency is higher than 1000 colonies per sgRNA construct in the library, harvest cells and purify the library

- 6) Depending on your pellet weight use multiple Maxiprep columns (Sigma or Zymogen), a Megaprep (Qiagen or Sigma), or Gigaprep (Qiagen, Sigma, Zymogen).
- 7) Expect a yield ~2mg.

Illumina sequencing

We strongly recommend deep sequencing the amplified libraries before use. Follow protocol "Illumina sequencing for amplified library" for optimized protocol on PCR-amplifying the sgRNA region and purification.